# MediCult IVM<sup>®</sup> System

Product No.:

8221

Languages:

## Customer Service:

E-mail: [customer.service@origio.com](mailto:customer.service@origio.com)  
Tel: +45 46 79 02 02 · Fax: +45 46 79 03 02

**ORIGIO a/s**  
Knardrupvej 2, DK-2760 Måløv, Denmark  
[www.origio.com](http://www.origio.com)  
Tel.: +45 46 79 02 00 Fax: +45 46 79 03 00

a CooperSurgical Company

## Ingredientia

**8201 LAG Medium – Vial 1**  
Albumini humani solutio; Aqua; CaCl<sub>2</sub>; Gentamicinum; Glucosum; KCl; MgSO<sub>4</sub>; NaCl; NaH<sub>2</sub>PO<sub>4</sub>; NaHCO<sub>3</sub>; Pyruvas; Phenolsulfophthaleinum; SSR

## 8220 IVM<sup>®</sup> Medium – Vial 2

Acetas; Acida aminica: Ala, Arg, Asp, Cys, Gln, Glu, Gly, His, Ile, Leu, Lys, Met, Phe, Pro, Ser, Thr, Trp, Tyr, Val, Cystinum, Glutathionum, Hydroxyprolinum; Adeninum; Adenosini phosphas; Adenosini triphosphas; Aqua; CaCl<sub>2</sub>; Cholesterolum; Desoxy-D-ribosum; D-ribosum; Fe(NO<sub>3</sub>)<sub>3</sub>; Gentamicinum; Glucosum; Guaninum; Hypoxanthinum; KCl; MgSO<sub>4</sub>; NaCl; NaH<sub>2</sub>PO<sub>4</sub>; NaHCO<sub>3</sub>; Phenolsulfophthaleinum; Polysorbatum 80; Pyruvas; Thyminium; Uracilium; Vitamini: A, B<sub>1</sub>, B<sub>2</sub>, B<sub>3</sub>, B<sub>5</sub>, B<sub>6</sub>, C, D<sub>2</sub>, E, Acidum 4-aminobenzoicum, Acidum nicotinicum, Cholinum, Inositolum, Menadionum, Nicotinamidum, Pyridoxalum, Pyridoxinum; Xanthinum

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en - english
--------------

For pre-incubation and maturing of immature oocytes.

This product is for ART treatment, whether the cause of infertility is male or female. The product should only be used by professionals trained in ART treatment.

**Package** 8201 LAG Medium (Vial 1) 8220 IVM® Medium (Vial 2)

**Pack size** 1 x Vial 1 (10 ml) 3 x Vial 2 (10 ml)

**Contains** Human albumin solution (HAS) (only Vial 1) Recombinant human insulin (only Vial 1) Gentamicin sulphate 10 µg/ml

**Quality control testing** Sterility tested (Ph.Eur., USP) Osmolality tested (Ph.Eur., USP) pH tested (Ph.Eur., USP) Endotoxin tested ≤ 0.1 EU/ml (Ph.Eur., USP) HSA analysis (Ph.Eur., USP) Mouse Embryo Assay (MEA) tested **Note:** The results of each batch are stated on a Certificate of Analysis, which is available at www.origio.com.

**Storage instructions and stability**
The products are aseptically processed and supplied sterile. Store in original container at 2-8°C, protected from light. Do not freeze. Discard excess (unused) media following warming. IVM® Medium is provided in vials intended for single use. LAG Medium is to be used within 7 days after opening. When stored as directed by the manufacturer the product is stable until the expiry date shown on the label.

**Precautions and warnings**
Do not use the product if:
1. Product packaging appears damaged or if the seal is broken.
2. Expiry date has been exceeded.
3. The product becomes discoloured, cloudy, or shows any evidence of microbial contamination.

**Caution:** The product is for single use and is not to be reused due to risk of contamination.

**Caution:** All blood products should be treated as potentially infectious. The source material used to manufacture this product was tested and found non-reactive for HBsAg and negative for Anti-HIV-1/-2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.

**Note:** Please note the need for traceability of this product. In addition, national legal requirements for this field may exist in your

country.

**Note:** Only to be used in combination with other devices intended for the particular purpose.

**Note:** Dispose of the device in accordance with local regulations for disposal of medical devices.

**Instructions for use**

- Pre-equilibrate 3 ml LAG Medium (Vial 1) and 10 ml of IVM® Medium (Vial 2) in a 5-6% CO<sub>2</sub> environment at 37°C for a minimum of 12 hours.
- After oocyte pick up, store the immature oocytes in LAG Medium in a 5-6% CO<sub>2</sub> environment at 37°C for 2-3 hours prior to transfer to the final IVM® Maturation Medium.
- Preparation of final IVM® Maturation Medium: 9 ml IVM® Medium (Vial 2) 1 ml patient's own serum 10 µl human chorionic gonadotropin (hCG) solution (100 mIU/ml) 100 µl human follicle stimulating hormone (FSH) solution (75 mIU/ml).
- Transfer the oocytes to the final IVM® Maturation Medium and incubate in a 5-6% CO<sub>2</sub> environment at 37°C for 28-32 hours.

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