English

Human Serum Albumin (HSA)

For laboratory procedures only; other uses must be qualified by the end user.

Product Description	REF Number	Unit Size
Human Serum Albumin 100 mg/mL	ART-3001	12 x 5 mL
Human Serum Albumin 100 mg/mL	ART-3003	100 mL

INTENDED USE

Human Serum Albumin is intended for use in assisted reproductive procedures which include gamete and embryo manipulations. This product is specifically intended for use as a protein supplement for SAGE™ reproductive culture media Quinn's Advantage™ Fertilization (HTF) Media, Ref. No. ART-1020 and ART-1021.

PRODUCT DESCRIPTION

Protein in the form of albumin, which makes up the bulk of the protein present in blood serum and reproductive tract fluid, is added to tissue culture media because it is thought to maintain the stability of cell membranes and chelate trace amounts of toxic components present in the culture water, media components, and containers and culture dishes. Protein, in the form of patient's serum or albumin, has been used extensively in media for IVF, GIFT, ICSI, embryo culture, embryo transfer, cryo-preservation, and sperm washing for IUI.

PRECAUTIONS AND WARNINGS

Do not use medium that shows evidence of particulate matter, or cloudiness.

The Human Serum Albumin (ART-3001 and ART-3003) contains 100 mg/mL of albumin in normal saline solution.

Caution: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when testing for antibodies to HIV-1/ HIV-2, HCV and non-reactive for HBsAg, HCV RNA and HIV-1 RNA. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents. Donors of the source material have been screened for Creutzfeldt-Jacobs disease (CJD). Based on effective donor screening and product manufacturing processes, it carries an extremely remote risk for transmission of viral diseases. A theoretical risk for transmission of CJD is also considered extremely remote. No cases of transmission of viral diseases or CJD have ever been identified for albumin.

Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes.

Single use: To avoid problems with contamination, handle using aseptic techniques and discard any excess product that remains in the bottle or vial after procedure is completed.

Reproductive media products are intended for single use only. Re-use of reproductive media may result in using a product past its labeled expiration date or increase the risk of microbial contamination in a subsequent procedure if the practitioner fails to utilize adequate aseptic techniques.

Use of expired or microbial contaminated product may result in suboptimal conditions to promote fertilization and/or embryo quality during in-vitro culture. These conditions may result in failure of the embryo to develop properly or to implant, potentially leading to a failed assisted reproductive procedure.

Note: Embryo is considered a general term. More precisely, SAGE™ considers the period of time initiating when a single diploid cell results from the fusion of male and female genome resulting in zygote formation with subsequent development from repeated mitotic divisions forming a solid mass or morula (typically day 4-5) and after which a fluid-filled cavity develops resulting in blastocyst formation (typically day 5-6) ending with embryo implantation that begins the end of the first week and is completed by the end of the second week post conception.

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

QUALITY ASSURANCE

Each lot of HSA is tested for pH (7.4 ± 0.2), osmolality (280 ± 10 mOsm/kg water), sterility (no detectable contamination).

One-cell MEA tested and passed with 80% or greater blastocyst. USP Endotoxin tested and passed with <1 EU/mL. A Certificate of Analysis is available for this product.

DIRECTIONS FOR USE

Quinn's Advantage™ Protein Plus Fertilization Medium, ART-1520, does not require protein supplementation as it contains 3 mg/mL of human serum albumin. Quinn's Advantage™ Fertilization Medium, ART-1020/ART-1021, is a protein-free medium that requires protein supplementation at the time of use. SAGE™ Human Serum Albumin 100 mg/mL Solution in Normal Saline is intended to provide the required protein supplementation to Quinn's Advantage™ Fertil+ization Medium, ART-1020/ART-1021.

When supplemented according to Table 1, Quinn's Advantage Fertilization Medium, ART-1020/ART-1021 will contain a final protein concentration of 3.0 mg/mL which is substantially equivalent to Quinn's Advantage™ Protein Plus Fertilization Medium, ART-1520.

TABLE 1

SUPPLEMENTATION OF QUINN'S ADVANTAGE™ FERTILIZATION MEDIUM, ART-1020/ART-1021		
Medium (mL)	Human Serum Albumin 100 mg/mL Solution (mL)	Final Volume (mL)
9.70	0.30	10.0
19.40	0.60	20.0
29.10	0.90	30.0
38.80	1.20	40.0
48.50	1.50	50.0
58.20	1.80	60.0
67.90	2.10	70.0
77.60	2.40	80.0
87.30	2.70	90.0
97.00	3.00	100.0

Refer to Quinn's Advantage™ Fertilization Medium, ART-1020/ART-1021 for further instructions for use.

NOT INTENDED FOR PARENTERAL USE IN HUMANS OR ANIMALS.

Information on specific aspects of IVF, embryo culture, and cryopreservation is available in our Product Catalog.

STORAGE INSTRUCTIONS AND STABILITY

Store unopened containers refrigerated at 2 °C to 8 °C. Warm to ambient or incubator (37 °C) temperature prior to use. Do not freeze or expose to temperatures greater than 39 °C. The product is stable until the expiration date shown on the label.

- Remove desired volume of product using aseptic procedures.
- Once removed, do not return any volume of product to the original container.
- Do not use if the product becomes discolored, cloudy, turbid, or shows any evidence of microbial contamination.

RELATED PRODUCTS

ART-1020/ART-1021 Quinn's Advantage™ Fertilization (HTF) Medium

SAGE™ In Vitro Fertilization™ has a full line of products for the Reproductive Medicine Specialist. Please call or write for specific information or to receive a copy of our current catalog. For technical questions, or to reach our Customer Service Department, call the SAGE™ Support Line.

Call the SAGE™ SUPPORT LINE: In the U.S.: (800) 243-2974 International: (203) 601-9818

EXPLANATION OF SYMBOLS

- REF** Catalog Number
- LOT** Batch Number
- USE BY** Use By (year, month, day)
- DO NOT REUSE** Do Not Reuse
- TEMPERATURE LIMITATION** Temperature Limitation
- STERILE A** Aseptic Technique Sterilization Membrane Filtered (SAL 10⁻³)
- ATTENTION:** See instructions for use.
- EC REP** Authorized Representative in the European Community.
- CE 0086** Product conforms to the Medical Device Directive 93/42/EEC
- MANUFACTURER** Manufacturer

RX ONLY U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

SAGE In Vitro Fertilization
a CooperSurgical Company
SAGE™ In Vitro Fertilization, Inc.
a CooperSurgical Company
95 Corporate Drive
Trumbull, CT 06611 USA
EC REP
ORIGIO a/s
Knardrupvej 2
DK-2760 Måløv
Denmark
www.origio.com

Customer Service:
E-mail: customer.service@origio.com
Tel: +45 46 79 02 02 Fax: +45 46 79 03 02

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