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1. Wilcox AJ, Baird DD, Weinberg CR. Time of implantation of the conceptus and loss of pregnancy. *N Engl J Med.* 1999;340:1796–1799
2. Garrido-Gómez T, et al. Profiling the gene signature of endometrial receptivity: clinical results. *Fertil Steril.* 2013;99:1078–1085
3. Teh WT, McBain J, Rogers P. What is the contribution of embryo-endometrial asynchrony to implantation failure? *J Assist Reprod Genet.* 2016;33:1419–1430
4. Ruiz-Alonso M, et al. The endometrial receptivity array for diagnosis and personalized embryo transfer as a treatment for patients with repeated implantation failure. *Fertil Steril.* 2013;100:818–824
5. Noyes RW, Hertig AT, Rock J. Dating the endometrial biopsy. *Fertil Steril.* 1950;1:561–564
6. Coutifaris C, et al. Histological dating of timed endometrial biopsy tissue is not related to fertility status. *Fertil Steril.* 2004;82:1264–1272
7. Ruiz-Alonso M, et al. What a difference two days make: “personalized” embryo transfer (pET) paradigm: a case report and pilot study. *Hum Reprod.* 2014;29:1244–1247
8. Petek M. qPCR, Microarrays or RNA-sequencing: When to choose one over the other? Available at: <https://biosistemika.com/blog/qpcr-microarrays-rna-sequencing-choose-one/> Last accessed June 2019
9. DAVID Bioinformatics Resource. Available at: <https://david.ncifcrf.gov/> Last accessed June 2019
10. Enciso M, et al. Development of a new comprehensive and reliable endometrial receptivity map (ER Map/ER Grade) based on RT-qPCR gene expression analysis. *Hum Reprod.* 2018;33:220–228
11. The Gene Ontology Resource. Available at: <http://geneontology.org/> Last accessed June 2019

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©2024 CooperSurgical Japan K.K. All Rights Reserved.

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Tel. 045-319-6826 Fax. 045-319-6581

E-mail: contact.genomics-jp@coopersurgical.com

Web: <https://coopersurgicalfertility-jp.com/>

CooperSurgical[®]
Fertility Solutions