




IUI BY ORIGIO®

# ORIGIO IUI Catheters

## Insemination Catheters



Single eyelet



Side eyelets



End eyelet

CE  
0543

origio

a CooperSurgical Fertility Company

## Versatile range to match the patient's anatomy and your preference

- Ergonomically designed handle for greater control
- Available in 3 different lengths
- Single and side eyelet options
- Gamma sterilized



Technical Information	ORIGIO IUI Catheters				
Catheter type	175mm	175mm Stiff	175mm Soft	140mm	60mm
Sheath length	135mm	135mm	135mm	100cm	N/A
Inner catheter (ID)	1.5mm	1.5mm	0.7mm	0.7mm	0.7mm
Inner catheter (OD)	1.9mm	1.9mm	2.0mm	2.0mm	2.0mm
Outer sheath (OD)	2.52mm	2.52mm	2.52mm	2.52mm	N/A
Volume of catheter not including hub	0.309mm	0.309mm	0.067mm	0.054mm	0.023mm
Volume per cm of catheter	0.018mm	0.018mm	0.004mm	0.004mm	0.004mm

### Product Compatibility and Order Codes:

	Artificial Insemination Catheter				
	INCA060 Quantity: 25 60mm	INCA140 Quantity: 25 140mm	INCA175 Quantity: 25 175mm	INCA175-O Quantity: 25 175mm	INCA175-OS Quantity: 25 175mm
Compatible with:					
2 side eyelets			✓		
End eyelet				✓	✓
Side eyelet	✓	✓			
Stiff				✓	
Soft					✓

✓ = Compatible

Legal Manufacturer: Smiths Medical International Ltd, Boundary Road, Hythe, Kent CT21 6JL, UK

**ORIGIO a/s**  
 Knardrupvej 2  
 2760 Måløv  
 Denmark  
 Tel: +45 46 79 02 00  
 www.origio.com



CooperSurgical Fertility Companies