



EMBRYO TRANSFER BY ORIGIO®

ORIGIO Spherical Tip

Supported Embryo Transfer Catheters



EchoGen™ technology
for visibility under
ultrasound



Spherical shaped
outer catheter



Hub designed
to ease handling
and grip

CE
0543

origio

a CooperSurgical Fertility Company

Our range of Spherical Tip catheters come individually sterile pouch packed for a 'one-step' or 'two-step' embryo transfer procedure.

- Spherical shape at the end of the outer catheter for greater control during catheter insertion
- Pre-formed angled tip 2cm from end of the outer catheter to aid insertion
- Unique hub for greater handling, grip, accuracy and control
- Soft formable stylet that does not compromise the performance of the outer sheath when in use
- Available in 24cm length
- Cervical stopper placed on the outer catheter to aid placement
- Available in our EchoGen range for visibility under ultrasound

Technical Information	ORIGIO Spherical Tip
Inner catheter (ID)	0.5mm
Inner catheter (OD)	1.00mm
Outer sheath (OD)	2.15mm
Volume of catheter not including hub	0.047cm ³ /cm
Volume per cm length of catheter	0.002cm ³ /cm
Stylet (OD)	1.10mm
Bulb size	2.6mm

Product Compatibility and Order Codes:

	EchoGen		Transem Classic	
	TSET24SASEG Quantity: 25 Soft 24cm Ultrasound	TSET24SASOEG Quantity: 25 Soft 24cm Ultrasound	TSET24SAS Quantity: 25 Soft 24cm Classic	TSET24SASO Quantity: 25 Soft 24cm Classic
Compatible with:				
Obturator included within pack		✓		✓
TSS1924 Quantity: 25. Transem Stylet for introducing the 24cm	✓	✓	✓	✓

✓ = Compatible

Legal Manufacturer: Smiths Medical International Ltd, Boundary Road, Hythe, Kent CT21 6JL, UK

ORIGIO a/s
Knardrupvej 2
2760 Måløv
Denmark
Tel: +45 46 79 02 00
www.origio.com



CooperSurgical Fertility Companies