



EQUIPMENT BY K-SYSTEMS®

K-Systems L200 Series

Class 2 sample and operator protection
(complies with EN12469)



a CooperSurgical Fertility Company

Features included as standard:

- Alu Heat technology, 5 heating zones – EDGE Enhancement
- Heated stainless steel surface with superior temperature control
- Silent Air Flow technology
- Clean working environment with HEPA filter
- Heated glass Ø90mm – LED light
- Ergonomic Danish design with Android touch controls

Options:

- Integrated RI Witness
- Several stand options including electric adjustment
- Monitor

Data	L234
Width (mm)	1246
Depth (mm)	735
Height (mm, standard stand)	2020
Power (kW)	0.9
Alu Heat technology	Standard
No. of heating zones	5

Heated Glass

Position	R
Size (mm)	Ø90
Separate heating zones	1
Light source	LED
Electrical outlets	1
Stereomicroscope preparation	Standard
Eye draught protection	Standard
Monitor (inches)	19 optional
Gassing/humidification	GH01 optional
120–240V	X

Filter

HEPA14	Standard
Flow sensor	Standard
RI Witness enabled	Optional
Support stand	Standard
Working height (mm, standard stand)	840
Keyboard shelf	Optional
PC shelf	Optional
Electrically adjustable height	Optional
Working height (mm, electrical stand)	790–1110
Packing included	Standard

Learn more about our equipment at www.origio.com