



EQUIPMENT BY K-SYSTEMS®

G185 Long Term Flat Bed Incubator

From the pioneers in Ambient-Air long-term incubation



1000+
units sold
worldwide

SYSTEMS
K

a CooperSurgical Fertility Company

Features included as standard:

- Reliable and stable gas and temperature control
 - HEPA and VOC filtered gas supply
- Built-in gas mixer – recirculated gas means low running costs
- Alu Heat technology heats the base plate
- Extensive data logging and monitoring capabilities
- Innovative Danish design

Data	G185 Standard	G185 ST
Width (mm)	850	850
Depth (mm)	558	558
Height (mm, standard stand)	152	152
Weight (kg)	38	38
Power (kW)	0.24	0.24
Gas consumption (l/h):		
CO ₂	0.6	0.6
N ₂	6	6
Recovery time (min):		
CO ₂	< 2	< 2
O ₂	< 4	< 4
Individual chambers	10	10
Equilibrium chamber	1	1
Temperature controlled	Standard	Standard
Gas controlled	No	No
Alu Heat technology	Standard	Standard
No. of heating zones	11	11
No. of temperature sensors	10	10
Dish inserts (min 10 per unit):		
NUNC	Optional	Optional
Falcon	Optional	Optional
Vitrolife	Optional	Optional
Continuous pH measurements		Optional
External CO ₂ measurements		Optional
External temperature measurements		Optional
SMS-enabled alerts	Standard	Standard
120–240V	X	X
Packing included	Standard	Standard

Options:

- Continuous pH measurements
- Independent external monitoring of temperature and CO₂
- Stacking system for up to 30 patients on same footprint



Stacking System	
Width (mm)	1100
Depth (mm)	830
Height (mm, standard stand)	1620
Weight (kg)	100

- Space for 3 incubators
- Space for PC
- Built-in power outlets
- Sliding shelves with soft-close feature
- Integrated ventilation system
- Drag chain to manage tubing and cables

The G185 is a CE marked device and is in conformity with the essential requirements of the Medical Devices Directive 93/42/EEC as amended by Directive 2007/47/EC.

Learn more about our equipment at www.origio.com

Legal Manufacturer: CooperSurgical Inc, 95 Corporate Drive, Trumbull, CT 06611, USA

ORIGIO a/s
 Knardrupvej 2
 2760 Måløv
 Denmark
 Tel: +45 46 79 02 00
www.origio.com



CooperSurgical Fertility Companies