



EQUIPMENT BY K-SYSTEMS®

# G185 Long Term Flat Bed Incubator

From the pioneers in Ambient-Air long-term incubation



1000+  
units sold  
worldwide



a CooperSurgical Fertility Company

## Features included as standard:

- Reliable and stable gas and temperature control
  - HEPA and VOC filtered gas supply
- Built-in gas mixer – recirculated gas means low running costs
- Alu Heat technology heats the base plate
- Extensive data logging and monitoring capabilities
- Innovative Danish design

Data	G185 Standard
Width (mm)	850
Depth (mm)	558
Height (mm, standard stand)	152
Weight (kg)	38
Power (kW)	0.24
Gas consumption (l/h):	
CO <sub>2</sub>	0.6
N <sub>2</sub>	6
Recovery time (min):	
CO <sub>2</sub>	< 2
O <sub>2</sub>	< 4
Individual chambers	10
Equilibrium chamber	1
Temperature controlled	Standard
Gas controlled	No
Alu Heat technology	Standard
No. of heating zones	11
No. of temperature sensors	10
Dish inserts (min 10 per unit):	
NUNC	Optional
Falcon	Optional
Vitrolife	Optional
Continuous pH measurements	
External CO <sub>2</sub> measurements	
External temperature measurements	
SMS-enabled alerts	Standard
120–240V	X
Packing included	Standard

## Options:

- Independent external monitoring of temperature and CO<sub>2</sub>
- Stacking system for up to 30 patients on same footprint



Stacking System	
Width (mm)	1100
Depth (mm)	830
Height (mm, standard stand)	1620
Weight (kg)	100

- Space for 3 incubators
- Space for PC
- Built-in power outlets
- Sliding shelves with soft-close feature
- Integrated ventilation system
- Drag chain to manage tubing and cables

## Learn more about our equipment at [www.origio.com](http://www.origio.com)

Legal Manufacturer: CooperSurgical Inc, 95 Corporate Drive, Trumbull, CT 06611, USA

**ORIGIO Inc**  
**CooperSurgical**  
 95 Corporate Drive  
 Trumbull, CT 06611  
 USA  
 Tel: 1-800-648-1151  
[www.origio.com](http://www.origio.com)



CooperSurgical Fertility Companies