



# ORIGIO Micropipettes

- Micropipettes since 1992
- No compromise on quality
- Customer-driven product development



## ICSI Micropipettes

Type	Reference no.	Available angles (XX in ref. no.)	Distance tip to bend	Bevel	Inner diameter	Other
ICSI SlimLine	MIC-SLM-XX	0, 20, 25, 30, 35, 45	0.5 mm	35-40°	4.3-4.9 µm	Parallel Taper
ICSI SlimLine	MIC-SLM-XX-B1.0	25, 30, 35	1 mm	35-40°	4.3-4.9 µm	Parallel Taper
ICSI	MIC-35-XX	0, 20, 25, 30, 35, 45	0.5 mm	30-35°	5.0-5.7 µm	Standard Taper
ICSI	MIC-35-XX-B1.0	20, 30, 35	1 mm	30-35°	5.0-5.7 µm	Standard Taper
ICSI	MIC-SI-XX	0, 20, 25, 30, 35, 45	0.5 mm	35-40°	5.0-5.7 µm	Standard Taper
ICSI	MIC-SI-XX-B1.0	20, 30, 35	1 mm	35-40°	5.0-5.7 µm	Standard Taper
ICSI	MIC-50-XX	0, 20, 25, 30, 35, 45	0.5 mm	50°	5.0-5.7 µm	Standard Taper
ICSI	MIC-50-XX-B1.0	25, 30, 35	1 mm	50°	5.0-5.7 µm	Standard Taper
Large ICSI	MIC-8-XX	20, 25, 30, 35, 40	0.5 mm	50°	6.2-7.0 µm	Standard Taper
Large ICSI	MIC-9-XX	20, 25, 30, 35	0.5 mm	40-45°	6.6-7.9 µm	Standard Taper

## Holding Micropipettes

Type	Reference no.	Available angles (XX in ref. no.)	Distance tip to bend	Inner diameter	Outer diameter
Small Holding	MPH-SM-XX	0, 20, 25, 30, 35, 45	0.5 mm	15-20 µm	65-95 µm
Medium Holding	MPH-MED-XX	0, 20, 25, 30, 35, 45	0.5 mm	15-20 µm	95-120 µm
Large Holding	MPH-LG-XX	0, 20, 25, 30, 35	0.5 mm	25-30 µm	120-150 µm
Extra Large Holding	MPH-XLG-XX	0, 25, 35	0.5 mm	25-30 µm	150-180 µm

## Biopsy Micropipettes

Type	Reference no.	Available angles (XX in ref. no.)	Distance tip to bend	Bevel	Inner diameter	Other
Biopsy	MPB-BS-XX	20, 30, 35	0.5 mm	45°	13-15 µm	Beveled and Spiked
Biopsy	MPB-BP-XX	30, 35	0.5 mm	45°	13-15 µm	Beveled and Polished
Biopsy	MPB-FP-XX	0, 20, 30, 35	0.5 mm	-	13-15 µm	Flat and Polished
Biopsy	MBB-FP-XXS-XX	25, 30, 35	0.5 mm	-	18-22 µm	Flat and Polished
Biopsy	MBB-FP-XS-XX	25, 30, 35	0.5 mm	-	23-27 µm	Flat and Polished
Biopsy	MBB-FP-SM-XX	0, 20, 25, 30, 35, 45	0.5 mm	-	28-32 µm	Flat and Polished
Biopsy	MBB-BP-SM-XX	20, 25, 30, 35, 45	0.5 mm	45°	28-32 µm	Beveled and Polished
Biopsy	MBB-FP-M-XX	0, 20, 25, 30, 35, 45	0.5 mm	-	33-37 µm	Flat and Polished
Biopsy	MBB-BP-M-XX	0, 20, 25, 30, 35	0.5 mm	45°	33-38 µm	Beveled and Polished
Biopsy	MBB-FP-L-XX	0, 20, 25, 30, 35	0.5 mm	-	38-42 µm	Flat and Polished
Biopsy	MBB-BP-L-XX	0, 30, 35	0.5 mm	45°	38-43 µm	Beveled and Polished

## Assisted Hatching Micropipettes

Type	Reference no.	Available angles (XX in ref. no.)	Distance tip to bend	Inner diameter
Assisted Hatching	MAH-SM-XX	0, 20, 25, 30, 35	0.5 mm	8-10 µm
Assisted Hatching	MAH-LG-XX	20, 25, 30, 35, 45	0.5 mm	10-12 µm

## Zona Dissection Micropipettes

Type	Reference no.	Available angles (XX in ref. no.)	Distance tip to bend
Zona Dissection	PZD-XX	0, 20, 25, 30, 35, 45	0.5 mm

ORIGIO products:



0086

ORIGIO a/s

Knardrupvej 2

2760 Måløv

Denmark

Tel: +45 46 79 02 00

www.origio.com



CooperSurgical Companies