



# RI Micropipettes

- High quality micropipettes
- Incorporating TwistaPak™ easy to use packaging



## ICSI Micropipettes

Type	Reference no.	Available angles (XX in ref. no.)	Distance tip to bend	Bevel	Inner diameter
ICSI Spiked Injection	7-71-T4SIXX/20	30, 35	0.75 mm	30°	4 µm
ICSI Spiked Injection	7-71-T5SIXX/20	30, 35	0.75 mm	30°	5 µm
ICSI Spiked Injection	7-71-T4SIXXL/20	30, 35	1 mm	30°	4 µm
ICSI Spiked Injection	7-71-T5SIXXL/20	30, 35	1 mm	30°	5 µm
ICSI Non-Spiked Injection	7-71-T5NIXX/20	30, 35	0.75 mm	18-20°	5 µm
ICSI Non-Spiked Injection	7-71-T5NIXXL/20	30, 35	1 mm	18-20°	5 µm

## Holding Micropipettes

Type	Reference no.	Available angles (XX in ref. no.)	Distance tip to bend	Inner diameter	Outer diameter
Holding	7-71-IHXX/20	0, 20, 25, 30, 35	0.75 mm	15 µm	90-110 µm
Large Holding	7-71-LHXX/20	20, 25, 30, 35	0.75 mm	25 µm	120-130 µm
Large Holding	7-71-32HXX/20	30, 35	0.75 mm	32 µm	170-180 µm

## Biopsy Micropipettes

Type	Reference no.	Available angles (XX in ref. no.)	Distance tip to bend	Bevel	Inner diameter	Other
Biopsy	7-71-14SBXX/20	35	0.75 mm	40°	14 µm	Beveled and Spiked
Biopsy	7-71-14NBXX/20	35	0.75 mm	40°	14 µm	Beveled
Biopsy	7-71-14FBXX/20	30, 35	0.75 mm	-	14 µm	Flat
Biopsy	7-71-20FBXX/20	30, 35	0.75 mm	-	20 µm	Flat
Biopsy	7-71-25FBXX/20	30, 35	0.75 mm	-	25 µm	Flat
Biopsy	7-71-30FBXX/20	30, 35	0.75 mm	-	30 µm	Flat
Biopsy	7-71-35FBXX/20	30, 35	0.75 mm	-	35 µm	Flat
Biopsy	7-71-40FBXX/20	30	0.75 mm	-	40 µm	Flat

## Assisted Hatching Micropipettes

Type	Reference no.	Available angles (XX in ref. no.)	Distance tip to bend	Inner diameter
Assisted Hatching	7-71-AHXX/20	30	0.75 mm	10 µm

## Zona Dissection Micropipettes

Type	Reference no.	Available angles (XX in ref. no.)	Distance tip to bend
Zona Dissection	7-71-MDXX/20	30, 35	0.75 mm

Research Instruments Micropipettes are CE 0120

### Research Instruments Ltd.

Bickland Industrial Park,  
Falmouth, Cornwall, UK  
TR11 2LF

Tel: +44 (0) 1326 372 753

www.research-instruments.com



CooperSurgical Companies