


Stat III™ Andrology Stains

Three Step Stain For Differentiation of Sperm in Semen Analysis.



- Results in 30 seconds
- Differentiates immature sperm and white blood cells
- "Strict" sperm morphology
- Excellent resolution of head, neck, and tailpiece
- Long term archiving capabilities



Spermatids (Immature Sperm Cells) can be differentiated from polymorphonuclear white blood cells (WBC) by using the following criteria: WBC contain connecting nuclear bridges and have a granular cytoplasm; Immature Sperm Cells are more spherical in shape.

For ordering information call your ORIGIO Representative at
1-800-648-1151

Stat III™ Andrology Stains

DESCRIPTION		ITEM #
STAT III™ Andrology Stain Stat-Pak	(3 x 90 ml)	MID-85303
STAT III™ Andrology Stain Set	(3 x 500 ml)	MID-85316
STAT III™ Fixative	(1 x 500 ml)	MID-85317
STAT III™ Solution I	(1 x 500 ml)	MID-85318
STAT III™ Solution II	(1 x 500 ml)	MID-85319
STAT III™ Fixative	(1 x Gallon)	MID-85328
STAT III™ Solution I	(1 x Gallon)	MID-85329
STAT III™ Solution II	(1 x Gallon)	MID- 85330



For ordering information call your ORIGIO Representative at
1-800-648-1151