







# The Stripper<sup>®</sup> and EZ-Range

- The EZ-Range is an innovative range of options for denudation and handling
- The Stripper; Your trusted brand for denudation handling



PIPETTORS	Type	Description	Reference no.
	EZ-Grip®	EZ-Grip micropipettor	7-72-2800
	EZ-Grip	EZ-Grip micropipettor with Comfort Rest	7-72-2802
	The Stripper	The Stripper micropipettor	MXL3-STR
	The Stripper	The Stripper micropipettor with Comfort Grip	MXL3-STR-CGR
	The Stripper BP	Bulb pipettor, black ergonomic	MXL3-STR-BP-DE
	The Stripper BP	Bulb pipettor, stainless steel slim	MXL3-STR-BP-SW
	The Stripper CC	Micropipettor for cumulus complex transfer	MXL3-STR-CC
	The Stripper PGD	Micropipettor for handling blastomeres	MXL3-STR-PGD
	The Stripper PGD	Micropipettor for handling blastomeres with Comfort Grip	MXL3-STR-PGD-CG

EZ-TIP®		Size of ID (µm)	Multipack (1x vial of 20)	Multipack (5x vial of 20)	Individual Pack (Box of 20)	Primary Use	Pipettor Options/Compatibility	
		EZ-Tip Vial	75	7-72-2075/1	7-72-2075/5	7-72-4075/20	Handling blastomeres for PGD /PGS or smaller species oocytes (eg mouse)  Oocyte stripping for ICSI	Stripper MXL3-STR; MXL3-STR-CGR
			125	7-72-2125/1	7-72-2125/5	7-72-4125/20		
			135	7-72-2135/1	7-72-2135/5	7-72-4135/20		
			145	7-72-2145/1	7-72-2145/5	7-72-4145/20		
			155	7-72-2155/1	7-72-2155/5	7-72-4155/20		
		EZ-Tip Individual	170	7-72-2170/1	7-72-2170/5	7-72-4170/20	Handling embryos	EZ-Grip with Comfort Rest 7-72-2800
			200	7-72-2200/1	7-72-2200/5	7-72-4200/20	Handling blastocysts and expanded blastocysts	EZ-Grip without Comfort Rest 7-72-2802
			290	7-72-2290/1	7-72-2290/5	7-72-4290/20	Handling larger cell masses and cumulus complexes	
			600	7-72-2600/1	7-72-2600/5	7-72-4600/20		

STRIPPER® TIPS		Size of ID (µm)	Multipack (1x vial of 20)	Individual Pack (Box of 20)	Primary Use	Pipettor Options/Compatibility	
		Stripper Tip Vial	50	MXL3-50	MXL3-IND-50	Handling blastomeres for PGD /PGS or smaller species oocytes (eg mouse)  Oocyte stripping for ICSI	Stripper PGD, MXL3-STR-PGD
			75	MXL3-75	MXL3-IND-75		
			100	MXL3-100	MXL3-IND-100		
			125	MXL3-125	MXL3-IND-125		
			135	MXL3-135	MXL3-IND-135		
		Stripper Tip Individual	150	MXL3-150	MXL3-IND-150	Handling embryos	EZ-Grip with Comfort Rest 7-72-2800
			175	MXL3-175	MXL3-IND-175	Handling blastocysts and expanded blastocysts	EZ-Grip without Comfort Rest 7-72-2802
			200	MXL3-200	MXL3-IND-200		
			275	MXL3-275	MXL3-IND-275	Handling larger cell masses and cumulus complexes	Stripper CC MXL3-STR-CC
600			MXL3-600	MXL3-IND-600			
		1000	MXL3-CC	MXL3-IND-CC			

ACCESSORIES	Type	Description	Reference no.
	EZ-Grip spares kit	Plunger, spacer and seal replacements	5-72-2000
	The Stripper O-ring repair kit	O-rings, spacer	MXL3-ORK
	The Stripper Plunger replacement	Plunger replacement kit	MXL3-PLG
	The Stripper CC O-ring repair kit	O-rings, spacer	MXL3-ORK-CC
	The Stripper PGD O-ring repair kit	O-rings, spacer	MXL3-ORK-PGD
	The Stripper PGD Plunger replacement	Plunger replacement kit	MXL3-PLG-PGD
	The Comfort Grip	Comfort Grip for The Stripper	MXL3-CGR
EZ-Stand	Stand for pipettors	6-72-2810	