



The Stripper[®]

Your trusted brand for denudation and handling



origio
'S

a CooperSurgical Company

Pipettors

| Type | Description | Reference no. |
|------------------|--|-----------------|
| The Stripper | The Stripper micropipettor | MXL3-STR |
| The Stripper | The Stripper micropipettor with comfort grip | MXL3-STR-CGR |
| The Stripper CC | Micropipettor for cumulus complex transfer | MXL3-STR-CC |
| The Stripper PGD | Micropipettor for handling blastomeres | MXL3-STR-PGD |
| The Stripper PGD | Micropipettor for handling blastomeres with comfort grip | MXL3-STR-PGD-CG |

Stripper Tips

| | Size of ID (µm) | Multipack (1x vial of 20) | Individual Pack (Box of 20) | Primary Use | Pipettor Options/Compatibility | |
|---|---|---------------------------|-----------------------------|--|---|--------------------------------|
|  Stripper Tip Vial | 50 | MXL3-50 | MXL3-IND-50 | Handling blastomeres for PGD/PGS | Stripper PGD, MXL3-STR-PGD | |
| | 75 | MXL3-75 | MXL3-IND-75 | Handling blastomeres for PGD/PGS or smaller species oocytes (e.g. mouse) | | |
| | 100 | MXL3-100 | MXL3-IND-100 | | Oocyte stripping for ICSI | |
| | 125 | MXL3-125 | MXL3-IND-125 | | | |
| | 135 | MXL3-135 | MXL3-IND-135 | | | |
| |  Stripper Tip Individual | 150 | MXL3-150 | MXL3-IND-150 | Handling embryos | StripperMXL3-STR; MXL3-STR-CGR |
| | | 175 | MXL3-175 | MXL3-IND-175 | | |
| | | 200 | MXL3-200 | MXL3-IND-200 | Handling blastocysts and expanded blastocysts | |
| | | 275 | MXL3-275 | MXL3-IND-275 | | |
| | | 600 | MXL3-600 | MXL3-IND-600 | Handling larger cell masses and cumulus complexes | |
| 1000 | | MXL3-CC | MXL3-IND-CC | | | |

Accessories

| Type | Description | Reference no. |
|--------------------------------------|-------------------------------|---------------|
| The Stripper O-ring repair kit | O-rings, spacer | MXL3-ORK |
| The Stripper Plunger replacement | Plunger replacement kit | MXL3-PLG |
| The Stripper CC O-ring repair kit | O-rings, spacer | MXL3-ORK-CC |
| The Stripper PGD O-ring repair kit | O-rings, spacer | MXL3-ORK-PGD |
| The Stripper PGD Plunger replacement | Plunger replacement kit | MXL3-PLG-PGD |
| The Comfort Grip | Comfort grip for The Stripper | MXL3-CGR |

