



The Original

Wallace V-Series ICSI Injection Pipette

For intracytoplasmic sperm injection



For more information about the Wallace V-tip and other micro-tools available, please contact your Wallace representative by visiting www.wallace-icsi.com

wallace[®]

a CooperSurgical Company

ICSI Injection Pipettes

Advantages



- Smooth control of sperm injection and cytoplasm suction
- Designed to be less damaging to the oocyte which may reduce cytoplasm draw back; this may increase the likelihood of better fertilisation rates and improve embryo development
- Multiple sperm can be loaded and injected

Characteristics



- Short taper
- Long parallel wall section
- Standard tip to elbow length (550µm)

Customer Testimonials

"We have always found them excellent quality and performed outstanding in every procedure."

Dr Levent Keskin-tepe, PhD,
HCLD, Executive Lab Director,
Sher Inst USA

"We have seen an improvement in fertilization rate and a decrease in damage rate."

C Franklin, Embryologist,
Nuffield, UK

"We are happy to report a reduction in lysis rate post-ICSI injection greater than 1% alongside continued pregnancy rates."

Dr Marta Jansa-Perez,
Head of Laboratories,
IVF Hammersmith

"Needle to needle consistency."

M Blayney,
Head of Science,
Bourn Hall, UK

Reference	Inner Diameter	Tip to elbow Lengths
WIC-45V – (angle)*	4.5µm	0.55mm**
WIC-50V – (angle)*	5.0µm	0.55mm**
WIC-55V – (angle)*	5.5µm	0.55mm**

*Angles available from 0-45 degrees

** 1.0mm Tip to elbow lengths available eg code reference WIC-45V-30L

ORIGIO a/s

Knardrupvej 2

2760 Måløv

Denmark

Tel: +45 46 79 02 00

www.origio.com



CooperSurgical Companies